
kasios[®]

TCH

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safe and efficient bone filling

The great number of studies and papers published on both HA (hydroxyapatite) and β TCP tricalcium phosphate show if need be the interest of such materials. For a long time these were the only calcium phosphates tested in vivo and implanted in humans. Their chemical composition close to that of bone favours biocompatibility and induces optimal osteointegration. The use of such interconnected porous biphasic ceramics (HA/TCP), as synthetic bone graft substitutes is one of the answers to **safe and efficient** repair of bone defects and for bone grafting.

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Composition

TCH[®] is a macroporous biphasic resorbable ceramic composed of 75% hydroxyapatite and 25% tricalcium phosphate.

Porosity - pore size




The TCH[®] totally interconnected porosity ranges from 60 to 80%. The mean pore diameter of 200-500 μ m is compatible with human bone-cell size.

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TCH[®] biocompatibility / biocompatibilité

Implanted in sheep condyles these ceramics have evidenced excellent integration in both cortical and cancellous bone.

1. Osteoconduction

After 6 weeks, a centripetal progressing neofomed bone can be observed enveloping the surface of the TCH[®] graft and having progressed through the interconnected macroporous structure to the centre of the graft.

2. Osteointegration

After 12 weeks, the TCH[®] graft is completely colonized by bone tissue covering the complete graft with no interposition of fibrous tissues.

3. Resorption - remodelling

After 12 weeks the macroporous structure of the graft is seen as a blurred outline and the resorption of the ceramic graft can be observed. Bone remodelling is a global phenomena; active osteoclasts are visible both at the ceramic surface and on the bone.

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A reference list of 406 clinical cases in different indications cases is available. After **over 10 years follow-up** these ceramics have resulted in entire satisfaction. Because of their synthetic origin no immunological or infectious reactions have ever been reported.

For system stabilization in the case of mechanical constraint the use of an osteosynthesis device is necessary. Clinical results are subordinated to good graft stability avoiding micro-movements, to optimal bone contact and to a healthy recipient site. Bone loss restoration has in every case been followed by rapid consolidation with perfect graft integration.

TCH® has been clinically successful in the following indications:

- for both epiphyseal and diaphyseal simple and complex fractures
- filling after removal of osteosynthesis materials and after benign synovioma curettage
- non-union or pseudarthrosis, arthrodesis and osteotomies
- prosthesis revision surgery
- spinal fusion

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1 an / 1 year

implantation / implantation

- Preferably implant TCH® in a cancellous bone site, after reaming to eliminate necrotic or inflammatory tissue.
- Impregnate TCH® with blood or even better with autologous marrow or bone.
- Fill to full capacity. Cover completely and hermetically.

As is the case for other porous ceramics, TCH® should not be subjected to mechanical stress in the absence of osteosynthesis.

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Shapes available:

- 2-3 mm irregular-shape porous granules. Particularly suitable for mixing with autologous bone and filling bony defects.
- Blocks, cylinders, wedges.

TCH® is gamma ray sterilized and packed in a double packaging. Do not re-sterilize. Single use product.

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benefits / avantages

- Synthetic (neither of human nor animal origin): no risk of infection or cross contamination.
- Bioresorbable, biocompatible.
- Excellent in-depth osteoconductive properties.

TCH® is a safe synthetic grafting material rapidly osteointegrated in recipient bone.

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P. Frayssinet, J.L. Trouillet, N. Rouquet, E. Azimus, A. Autefage. 1993. Effects of the chemical composition of calcium phosphates ceramics on their osseointegration. Orthopaedics International Edition. July/August/September Vol.1 No.4.

P. Frayssinet, J.L. Trouillet, N. Rouquet, E. Azimus, A. Autefage. 1993. Osseointegration of macroporous calcium phosphate ceramics having a different chemical composition. Biomaterials, Vol.14 No.6.

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